

Ultrasound Biopsy Kits  
Disposable Biopsy Guides  
Ultrasound Cover Kits  
Reusable biopsy guide

# Catalogue



## Products



### Ultrasound Biopsy Kits

- General Purpose Kit include: Needle Guide, Transducer Cover, Sterilized Gel & 2 Elastic Bands
- Innovative Guide Design, Totally Disposable, Less infection risk
- Full Range of Needle Adapters for different procedures/gauges
- Telescopically-folded , soft and durable PU material , latex-free transducer cover
- Single use Sterilized Gel reduces infection risk



### Disposable Biopsy Guides

- Multi Angle Guide, Free-Angle Guide & Endocavity Guide for different procedures
- Innovative Design, Totally Disposable, Less infection risk due to Single-Use
- Full Range of Needle Adapters for different procedures/gauges



### Ultrasound Cover Kits

- Cover kits include: Transducer Cover, Sterilized Gel & 2 Elastic Bands
- Single-Use, reduces the risk of cross-contamination
- Telescopically-folded , soft and durable PU material , latex-free Transducer cover



### Reusable biopsy guide

- Reusable stainless steel bracket and tubes/adapters
- Over 300 type of guides
- Can be customized according to special probe shape and procedure
- Medical stainless steel material 304 or 316L
- Wide range of needle gauge selection: 14G-23G , 6F-12F



## Ultrasound Biopsy Kits

- › Innovative Guide Design, Totally Disposable, Less infection risk
- › General Purpose Guide include: Needle Guide, Transducer Cover, Sterilized Gel & 2 Elastic Bands
- › Widely used needle guides offer accurate guidance for various applications including core tissue biopsy, tumor ablation, fluid aspiration, catheter placement, marker placement, regional anesthesia, vascular access



## Ultrasound Cover Kits Sterilized Packing

Single-Use, reduces the risk of cross-contamination

Telescopically-folded, soft and durable PU material,

Latex-free Transducer cover




## Usage Instructions

- 1) Place an appropriate amount of gel inside the cover or on the transducer lens. Lack of gel may result in a poor imaging.
- 2) Insert transducer into cover, and make sure to use sterile technique. Cover the transducer cable as far as possible.
- 3) Make sure there is no air between transducer lens and cover.
- 4) Pull the cover tightly over transducer. Tie the bands.
- 5) Inspect the cover to ensure there are no holes or tears.



Product	Model	Order No.	Description	Qty
Cover Kits	LPTCKG01	10500102	14*90cm PU cover,20ml Gel	24
	LPTCKG01	10500202	14*120cm PU cover,20ml Gel	24
	LPTCKG01	10500302	14*150cm PU cover,20ml Gel	24
	LPTCKG01	10500602	14*61cm PU cover,20ml Gel	24
	LPTCKG01	10500702	9*92cm PU cover,20ml Gel	24
Transducer Cover	LPCOVG	10400102	14*90cm PU cover only	50
	LPCOVG	10400202	14*120cm PU cover only	50
	LPCOVG	10400302	14*150cm PU cover only	50
	LPCOVG	10400602	14*61cm PU cover only	50
	LPCOVG	10400702	9*92cm PU cover only	50
Gel	LPGELSa20	10300102	20ml Sterilized Gel	100




## Multi-Angle Guide


Accurate Guidance for core or in-depth biopsy

- Multi-Angle Guide utilizes a two-part system consisting of a disposable guide and disposable needle adapters.
- Multi-Angle Guide provides a fixed path according to US machine software biopsy-line. Push the Angle-switch button to select the angle.
- Accepts 6F-10F & 14G-23G instruments

Multi-Angle Guide can be used in a variety of applications, such as biopsy, ablation, and drainage.



The guide



The adapters

## Step by Step Instructions:



**1. Install the cover**  
Apply gel onto the transducer, then put on the cover to reduce the risk of infection.



**2. Install the guide**  
Place the guide on the transducer with the needle adapter on the appropriate side, and put the lock in place to stabilize the guide.




**3. Select an adapter**  
Choose the appropriate adapter for your needle. Push the lateral buttons to switch the angle that matches your biopsy system's settings.



**4. Release the needle**  
During the procedure, you can release the needle by taking out the adapter.



**5. Dispose the guide**  
Pull the lock strip to release the guide, and take the guide off the transducer.




## Free-Angle Guide

Flexibly Guidance for Superficial Biopsy, Regional Anesthesia, Vascular access

- Free-Angle Guide utilizes a two-part system consisting a disposable guide and disposable needle adapters. Innovative all-in-one single-use.
- Free-Angle Guide provides infinite angle accurate in-plane, also allows an efficient release.
- Accepts 14G,16G,17G,18G,20G,22G instruments



The guide



The adapters


## Step by Step Instructions:



**1. Install the cover**  
Apply gel onto the transducer, then put on the cover to reduce the risk of infection.




**2. Install the guide**  
Place the guide on the transducer with the needle holder on the appropriate side, and put the lock in place to stabilize the guide.



**3. Select an adapter**  
Choose the appropriate adapter for your needle.



**4. Release the needle**  
During the procedure, you can release the needle from the adapter.




**5. Dispose the guide**  
Pull the lock strip to release the guide, and take the guide off the transducer


## Endocavity guide

Long Tunnel design for IVF and Prostate Biopsy  
Short Tunnel design for Prostate Biopsy

- Endocavity Guide provide a fixed path according to US machine software biopsy-line. The short tunnel guide is specially designed for prostate biopsy. The long tunnel guide is used for both prostate biopsy and IVF.
- Endocavity Guides are recommended for procedures such as: tissue biopsy, fluid aspiration, marker placement, placement of radioactive seeds.
- Accepts 16G-18G instruments



Endoscopic guide  
(long)




Endoscopic guide  
(short)

## Step by Step Instructions:



### 1. Install the cover

Apply gel onto the transducer, then put on the cover to reduce the risk of infection.



### 2. Install the guide

Place the guide on the transducer and align with the transducer attachment holes, put the lock in place to stabilize the guide.



### 3. Insert the needle

Insert the needle through the tunnel



### 4. Dispose the guide

Pull the lock strip to release the guide, and take the guide off the transducer



## Reusable biopsy guides

- Leapmed's reusable biopsy guides are non-sterilized.
- The reusable guide requires proper sterilization and disinfection before each use.
- Reusable guides are made of 304 or 316L stainless steel material.
- Standard General-use biopsy guide contains a stainless bracket and two tubes, 16G and 18G .  
Other tubes on demands from 14G-23G & 6F-12F.
- Endocavity Guide accepts 16G-18G instruments.







Stable-  
Adaptable-  
Disposable

Disposable Endocavity short guide



Compatible Probe	Order No.	Type	Model	Description	Tube	Qty/Box
GE 11L,11L-D,12L-RS,12L-SC,L4-12L-RS,L6-12-RS	10609702	Free-Angle	LPUBKG56	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-22G	12
GE 9L,9L-D,9L-RS	102046	Free-Angle	LPG9L-free	Guide only	14G-22G	10
	10604602	Free-Angle	LPUBKG19	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-22G	12
GE ML6-15-D, ML6-15-RS	102047	Free-Angle	LPGML615-free	Guide only	14G-22G	10
	10604702	Free-Angle	LPUBKG24	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-22G	12
GE 11L,11L-D,12L-RS,12L-SC,L4-12L-RS,L6-12-RS	102042	Guideline-Angle	LPG12L	Guide only	14G-23G,6F-10F	10
	10604202	Guideline-Angle	LPUBKG15	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
GE 3CRF-D	102043	Guideline-Angle	LPG3CRFD	Guide only	14G-23G,6F-10F	10
	10604302	Guideline-Angle	LPUBKG16	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
GE 4C,4C-A,4C-D,4C-RS,4C-SC	102011	Guideline-Angle	LPG4C	Guide only	14G-23G,6F-10F	10
	10601102	Guideline-Angle	LPUBKG17	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
GE 9L,9L-D,9L-RS	102018	Guideline-Angle	LPG9L	Guide only	14G-23G,6F-10F	10
	10601802	Guideline-Angle	LPUBKG18	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
GE C1-5-D, C1-5-RS	102004	Guideline-Angle	LPGC15	Guide only	14G-23G,6F-10F	8
	10600402	Guideline-Angle	LPUBKG25	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
GE C1-6-D,C1-6VN	102012	Guideline-Angle	LPGC16	Guide only	14G-23G,6F-10F	10
	10601202	Guideline-Angle	LPUBKG20	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
GE C2-6b-D	102053	Guideline-Angle	LPGC26bD	Guide only	14G-18G,6F-8.5F	10
	10605302	Guideline-Angle	LPUBKG21	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-18G,6F-8.5F	12
GE C2-7,C2-7-VN	10609302	Guideline-Angle	LPUBKG52	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
	102052	Guideline-Angle	LPGC29D	Guide only	14G-23G,6F-10F	10
GE C2-9-D	10605202	Guideline-Angle	LPUBKG22	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
	102010	Guideline-Angle	LPGML615	Guide only	14G-23G,6F-10F	10
GE ML6-15-D, ML6-15-RS	10601002	Guideline-Angle	LPUBKG23	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
	10208502	Long-Endo	LPERIC59-L	Guide only	16G-18G	24
GE RIC5-9-D, RIC6-12-D, RIC5-9W-RS	102051	Short-Endo	LPERIC59	Guide only	16G-18G	10
	102001	Long-Endo RU Model	LPEIC59	Guide only	16G-18G	10
GE E8C, E8C-RS, E8CS, IC5-9-D, IC5-9H	102066	Short-Endo TR5 Model	LPEIC59-short	Guide only	16G-18G	10

Stable-  
Adaptable-  
Disposable

Disposable Free angle guide



Compatible Probe	Order No.	Type	Model	Description	Tube	Qty/Box
Phillips L12-3(Affiniti70, EPI05, EPI07, CX50), L12-4(Affiniti 50, ClearVue350, 550, 650, 850)	10603502	Free-Angle	LPUBKG34	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-22G	12
	102035	Free-Angle	LPGL123CX	Guide only	14G-22G	10
Phillips L12-5	102022	Free-Angle	LPGL125	Guide only	14G-22G	10
	10602202	Free-Angle	LPUBKG36	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-22G	12
Phillips C5-1	102006	Guideline-Angle	LPGC51	Guide only	14G-23G,6F-10F	8
	10600602	Guideline-Angle	LPUBKG30	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Phillips C9-2	102041	Guideline-Angle	LPGC92	Guide only	14G-23G,6F-10F	8
	10604102	Guideline-Angle	LPUBKG31	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Phillips ClearVue C5-2	102036	Guideline-Angle	LPGC52CV	Guide only	14G-23G,6F-10F	10
	10603602	Guideline-Angle	LPUBKG32	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Phillips Ergo C5-2 (ENVISOR, HD5, HD6, HD7, HD7XE, HD11, HD11XE, HDI 4000, HDI 5000/5500cv, IE33,IU22)	10603402	Guideline-Angle	LPUBKG35	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
	102034	Guideline-Angle	LPGC52HD	Guide only	14G-23G,6F-10F	8
Phillips Ergo S4-2,S5-1,S5-2	10609602	Guideline-Angle	LPUBKG55	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Phillips L12-4( Affiniti 50, ClearVue 350,550,650,850)	102030	Guideline-Angle	LPGL124CV	Guide only	14G-23G,6F-10F	8
	10603002	Guideline-Angle	LPUBKG33	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Phillips BP10-5ec	10207602	Long-Endo	LPEE312A	Guide only	16G-18G	24
Phillips C8-4V,C10-3V	102003	Long-Endo	LPEC84V/103V	Guide only	16G-18G	10
	102062	Short-Endo	LPEC84V-short	Guide only	16G-18G	10
Phillips C9-4V(Affiniti,CX30,CX50,Sparq)	10207402	Long-Endo	LPEAFC94V	Guide only	16G-18G	24
	102029	Long-Endo	LPEC94V	Guide only	16G-18G	10
Phillips ClearVue C9-4V	102061	Short-Endo	LPEC94V-short	Guide only	16G-18G	10
	102027	Short-Endo	LPEC95	Guide only	16G-18G	10

Stable·  
Adaptable·  
Disposable

Disposable Angle-switch guide



Compatible Probe	Order No.	Type	Model	Description	Tube	Qty/Box
Esaote LA332,SL332	102031	Free-Angle	LPGLA332	Guide only	14G-22G	10
	10603102	Free-Angle	LPUBKG12	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-22G	12
Esaote LA522, LA523,LA532	102021	Free-Angle	LPGLA523	Guide only	14G-22G	10
	10602102	Free-Angle	LPUBKG13	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-22G	12
Esaote AC2541	102038	Guideline-Angle	LPGAC2541	Guide only	14G-23G,6F-10F	10
	10603802	Guideline-Angle	LPUBKG10	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Esaote AL2443, LA533,SL1543	102039	Guideline-Angle	LPGSL1543	Guide only	14G-23G,6F-10F	10
	10603902	Guideline-Angle	LPUBKG14	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Esaote CA541	102008	Guideline-Angle	LPGCA541	Guide only	14G-23G,6F-10F	10
	10600802	Guideline-Angle	LPUBKG11	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Esaote TRT33	102040	Short-Endo	LPETRT33	Guide only	18G	10
Esaote EC1123,SE3123,SE3133	102013	Short-Endo 0° model	LPEEC1123-0	Guide only	16G-18G	10
Esaote EC1123,EC123,SE3123,SE3133	102059	Short-Endo 3.8° model	LPEEC1123-3.8	Guide only	16G-18G	10
Siemens 9L4	102048	Free-Angle	LPG9L4-free	Guide only	14G-22G	10
	10604802	Free-Angle	LPUBKG41	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-22G	12
Siemens LA522E,LA523	102021	Free-Angle	LPGLA523	Guide only	14G-22G	10
	10602102	Free-Angle	LPUBKG13	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-22G	12
Siemens 18L6HD	102028	Guideline-Angle	LPG18L6	Guide only	14G-23G,6F-10F	10
	10602802	Guideline-Angle	LPUBKG37	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Siemens 4C1,6C1HD,CH5-2,CH4-1	102007	Guideline-Angle	LPG4C1/6C1	Guide only	14G-23G,6F-10F	8
	10600702	Guideline-Angle	LPUBKG38	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Siemens 4V1	102002	Guideline-Angle	LPG4V1	Guide only	14G-23G,6F-10F	10
	10600202	Guideline-Angle	LPUBKG39	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Siemens 9L4	102023	Guideline-Angle	LPG9L4	Guide only	14G-23G,6F-10F	10
	10602302	Guideline-Angle	LPUBKG40	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Siemens E3-12A, EC10-5w	10207602	Long-Endo	LPEE312A	Guide only	16G-18G	24
Siemens EC9-4(S1000,S2000,S3000,Antares)	102009	Long-Endo	LPEEC94S	Guide only	16G-18G	10
Siemens EC9-4,EV9-4(Sonoline G40,X150,X300)	102005	Long-Endo	LPEEC94X	Guide only	16G-18G	10
	102058	Short-Endo	LPEEC94X-short	Guide only	16G-18G	10
Siemens MC9-4,EC9-4w	102017	Long-Endo	LPEMC94	Guide only	16G-18G	10
	102013	Short-Endo 0° model	LPEEC1123-0	Guide only	16G-18G	10
Siemens EC1123	102059	Short-Endo 3.8° model	LPEEC1123-3.8	Guide only	16G-18G	10

Stable·  
Adaptable·  
Disposable


Disposable Endocavity guide



Compatible Probe	Order No.	Type	Model	Description	Tube	Qty/Box
Aloka UST-5415	102045	Guideline-Angle	LPGL441	Guide only	14G-23G,6F-10F	10
	10604502	Guideline-Angle	LPUBKG06	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Aloka UST-9130,9126,9123,979	102044	Guideline-Angle	LPG9130	Guide only	14G-23G,6F-10F	10
	10604402	Guideline-Angle	LPUBKG07	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Aloka UST-9133,C22K	102064	Guideline-Angle	LPG9133	Guide only	14G-18G,6F-10F	10
	10606402	Guideline-Angle	LPUBKG08	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-18G,6F-10F	12
Aloka UST-679	10207602	Long-Endo	LPEE312A	Guide only	16G-18G	24
Aloka UST-981-5,UST-984-5	102026	Long-Endo	LPE984	Guide only	16G-18G	10
Aloka UST-9124,9118	102037	Short-Endo	LPE9124	Guide only	16G-18G	8
Hitachi L74ML55	102049	Free-Angle	LPGL74M-free	Guide only	14G-22G	10
	10604902	Free-Angle	LPUBKG04	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-22G	12
Hitachi L75,L64	102069	Free-Angle	LPGL75	Guide only	14G-22G	10
	10606902	Free-Angle	LPUBKG05	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-22G	12
Hitachi B512	102065	Guideline-Angle	LPGB512	Guide only	14G-23G,6F-10F	10
	10606502	Guideline-Angle	LPUBKG01	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Hitachi C251	102056	Guideline-Angle	LPGC251	Guide only	14G-23G,6F-10F	10
	10605602	Guideline-Angle	LPUBKG09	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Hitachi C715,C25	102024	Guideline-Angle	LPGC715	Guide only	14G-23G,6F-10F	10
	10602402	Guideline-Angle	LPUBKG02	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Hitachi L441	102045	Guideline-Angle	LPGL441	Guide only	14G-23G,6F-10F	10
	10604502	Guideline-Angle	LPUBKG06	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Hitachi L74ML55	102025	Guideline-Angle	LPGL74M	Guide only	14G-23G,6F-10F	10
	10602502	Guideline-Angle	LPUBKG03	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Hitachi V53W/CC531	102014	Long-Endo	LPEV53W	Guide only	16G-18G	10
	102060	Short-Endo	LPEV53W-short	Guide only	16G-18G	10
Hitachi V73W,C41V1	102063	Long-Endo	LPEV73W	Guide only	16G-18G	10
	102055	Short-Endo	LPEV73W-short	Guide only	16G-18G	10

Stable·  
Adaptable·  
Disposable


Disposable Free angle guide



Compatible Probe	Order No.	Type	Model	Description	Tube	Qty/Box
Mindray L12-4S, L12-3E, 7L4A, 7L4	102050	Free-Angle	LPGL123E-free	Guide only	14G-22G	10
	10605002	Free-Angle	LPUBKG27	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-22G	12
Mindray C5-2, C5-2E, C5-2S	102019	Guideline-Angle	LPGC52E	Guide only	14G-23G,6F-10F	10
	10601902	Guideline-Angle	LPUBKG28	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Mindray L12-4S, L12-3E, 7L4A, 7L4	102020	Guideline-Angle	LPGL123E	Guide only	14G-23G,6F-10F	10
	10602002	Guideline-Angle	LPUBKG26	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Mindray SC5-1U	102068	Guideline-Angle	LPGSC51U	Guide only	14G-23G,6F-10F	10
	10606802	Guideline-Angle	LPUBKG29	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Mindray DE10-3U	10208502	Long-Endo	LPERIC59-L	Guide only	16G-18G	24
Mindray EV6 5MHz/R10, 6CV1, 6CV1s, 65EB10EA,V10-4,V10-4B, V10-4s, V11-3s,V11-3BE,V11-3EV11-3HU,V11-3Ws	102015	Long-Endo	LPEV104	Guide only	16G-18G	10
	102016	Short-Endo	LPEV104-short	Guide only	16G-18G	10
Mindray DE10-3U	102051	Short-Endo	LPERIC59	Guide only	16G-18G	10
Toshiba PVT-1005BT,805AT	102033	Guideline-Angle	LPG1005BT	Guide only	14G-23G,6F-10F	10
	10603302	Guideline-Angle	LPUBKG45	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Toshiba PVT-375BT	102032	Guideline-Angle	LPG375BT	Guide only	14G-23G,6F-10F	10
	10603202	Guideline-Angle	LPUBKG46	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Toshiba PVT-382BT	10609802	Guideline-Angle	LPUBKG57	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G,15G,16G,18G,20G,8.5F	12
Toshiba PVT-661VT,651VT	102009	Long-Endo	LPEEC94S	Guide only	16G-18G	10
Toshiba PVT-781VT,PVU-781VT,PVT-781VTE,PVU-781VTE	10207902	Long-Endo	LPE781VT	Guide only	16G-18G	24



Stable·  
Adaptable·  
Disposable



Disposable Free angle guide








Compatible Probe	Order No.	Type	Model	Description	Tube	Qty/Box
Alpinion E3-10	10207602	Long-Endo	LPEE312A	Guide only	16G-18G	24
BK 8823	10608902	Guideline-Angle	LPUBKG48	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
BK Medical 8819, E10C4	10207602	Long-Endo	LPEE312A	Guide only	16G-18G	24
AccuCARE Stepper	10207702	Short-Endo	LPEGRID	Guide only	17G-18G	24
Carestream E10C4	10207602	Long-Endo	LPEE312A	Guide only	16G-18G	24
Samsung CS1-4	102011	Guideline-Angle	LPG4C	Guide only	14G-23G,6F-10F	10
	10601102	Guideline-Angle	LPUBKG17	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-23G,6F-10F	12
Samsung Medison EVN4-9,ER4-9/10ED	102015	Long-Endo	LPEV104	Guide only	16G-18G	10
Samsung Medison NEV4-9,ESNER4-9ES	10207602	Long-Endo	LPEE312A	Guide only	16G-18G	24
Samsung VR5-9	102001	Long-Endo RU Model	LPEIC59	Guide only	16G-18G	10
	102066	Short-Endo TR5 Model	LPEIC59-short	Guide only	16G-18G	10
Samsung EV3-10B,V5-9	10208502	Long-Endo	LPERIC59-L	Guide only	16G-18G	24
	102051	Short-Endo	LPERIC59	Guide only	16G-18G	10
Samsung Medison EVN4-9,ER4-9/10ED	102016	Short-Endo	LPEV104-short	Guide only	16G-18G	10
SonoSite C60, C60e,C60x	102057	Free-Angle	LPGC60	Guide only	14G-22G	10
	10605702	Free-Angle	LPUBKG42	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-22G	12
SonoSite HFL38, HFL38x, L38, L38e and L38x	102070	Free-Angle	LPGHFL38	Guide only	14G-22G	10
	10607002	Free-Angle	LPUBKG43	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-22G	12
SonoSite HFL38xi, HFL38xp, L38xi,L38xp	10609402	Free-Angle	LPUBKG53	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-22G	12
Sonosite HFL50x, HFL50xp	10609002	Free-Angle	LPUBKG49	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	14G-22G	12
SonoSite L25, L25e, L25x,L25n,L25xp	102054	PICC	LPGL25N	Guide only	18G,21G (1.0/1.5/2/2.5/3,	10
	10605402	PICC	LPUBKG44	Biopsy Kits(Guide,20ml Gel,14*150cm PU Cover)	18G,21G (1.0/1.5/2/2.5/3,	12
SonoSite C8x	102027	Short-Endo	LPEC95	Guide only	16G-18G	10
Supersonic SE12-3,SEV12-3	10207802	Long-Endo	LPSE123	Guide only	16G-18G	24
Vinno G4-9E	10207802	Long-Endo	LPSE123	Guide only	16G-18G	24

Product list of Reusable guide

	
<p>C1-5-D, C1-5-RS 4C,4C-A,4C-D,4C-RS,CS1-4 C2-9-D ML6-15-D,ML6-15-RS E8C(LOGIQ C5) E8C,E8CS,E8C-RS,IC5-9-D 3Sc-RS S1-5-D (E9) 3CRF-D C2-6b-D 8C M5S 8L 3S E7C (GE LOGIQ C5) RIC5-9-D, RIC6-12-D, RIC5-9W-RS, V5-9 AB2-7 MTZ C2-7,C2-7-VN,C2-7-D RAB4-B-D 4C-RC (LOGIQ C5) 12L(LOGIQ3/5/7/9) C4-8 C1-6-D,C1-6-VN 9L,9L-D, 12L,12L-RS,12L-RC,L4-12L,L6-12-RS,11L,11L-D</p>	<p>C5-1 L12-5 L12-5 L12-3(HD Series) L12-3(CX50- EPIO5/7- Affinit50/70) L12-4(Sparq,CX30) L12-4(ClearVue350,550,650,850,Affinit50) C5-2(ClearVue350,550,650,850) C5-2 (HD Series- IU22) C8-4V,C10-3V,C9-3V,C9-5ec,C9-5 C9-4V(ClearVue) C10-4ec Ergo S4-2,S5-1,S5-2(CX,EPIQ) S5-1(IU22) X6-1 C9-4V(CX30,Affinit30/50/70) C6-2(Sparq,CX30,Affinit30/50/70) C8-5 L9-3 S4-1(HD Series,IU22,IU Elite) BP10-5EC L18-5 3D9-3V E6509</p>

	
<p>4C1,CH5-2,CH4-1 6C1HD EV-8C4 9L4 14L5 4V1 18L6HD EC9-4,EV9-4(X150- X300) MC9-4,EC9-4W 4V1c 6L3 8C3HD VF13-5 VF12-4</p>	<p>EUP-C715,C25 EUP-V53W,CC531 EUP-V73W,C41V1 EUP-B512 UST-9123,9130,9126 UST-5412,5413 UST-9124,9118 EUP-B715 EUP-L75,L64 UST-984-5,UST-981-5 UST-9147 UST-9120,C42K UST-9133 C532,C42 C514 B712 C251 L441,UST-5415 UST-678</p>

	
<p>CA541 LA522, LA523,LA532 LA332,SL332 Esaote LA332,SL332 CA431,CA421,CA621,CA631 EC1123,SE3123 Esaote EC1123,SE3123,EC123 TRT33 SL3323 PA240 CA123 SC3421 AL2443,SL1543,LA533 SE3133 SC3123 SP3630 AL2443, LA533,SL1543 AC2541</p>	<p>CA1-7A EV4-9ES L3-12A,LA3-14AD L5-13IS SC1-6 C2-6IC, EC65,EC4-9ED LA3-16A</p>

		
<p>C5-2,C5-2E,C5-2S 3C5A,3C5S 7L4A,7L4S,7L4,L12-3E,L12-4S V10-4,V11-3E,V11-3BE,6CV1,V11-3HU,6CV1S SC5-1U SC5-1E L14-5W C5-1U,C5-1S,C5-1E</p>	<p>HFL38x P21X L25, L25e, L25x, L25n, L25xp C60x, C60, C60e ICTx C60xi, C60xp</p>	<p>PVT-375BT PLT-805AT,1005BT PVT-661VT,PVT-651VT PVT-681MV PVT-770RT PST-30BT 781MV,781VT 704AT PVT-382BT PLT-1204BT PVT-625RT PVM-621VT PLT-712SBT,PVT-712BT PVM-375AT PLT-1204MV PLT-705BT</p>

		
<p>Chison D3C20L Chison D3C60L Chison D3C60L Chison D7C10L Chison D7L40L,Chison D7L40L Chison D7L40L Chison D6C15L Chison M8L50L mSonics L10-5 mSonics L10-5 Neusoft C2-5 Neusoft C4-8V Neusoft C6-2 Neusoft L12-4 Neusoft S4-2C Neusoft C8-4ec Neusoft C4-1 Vinnco G2-5C Vinnco G4-9E Vinnco G4-9E Vinnco S1-6P Vinnco X6-16L Vinnco F2-5C Vinnco G1-4P Vinnco S1-8C Vinnco F4-9E Sonoscape C344 Sonoscape L741 Sonoscape BCL10-5 Sonoscape 6V1 Sonoscape 6V3 Sonoscape 3C-A</p>	<p>Supersonic SE12-3 Supersonic SL10-2 Supersonic SC6-1 TCL C5-2/60S TCL L10-5/38S TCL EV9-4/S Telemed C3.5 Telemed C5-2R60HI Terason 12L5A SIUI L8L38C SIUI C3L60C</p>	

## LEAPMED HEADQUARTERS


Leapmed Healthcare Corporation  
1F,Block3,No.38,North Guandu Road  
CN-21510 Suzhou,Jiangsu,P.R.China  
Tel:(+86)512 66554498  
Fax(+86)512-68223292 68127433-8002

## Leapmed Europe

5 CHEMIN DE LABORIE  
31560 NAILLOUX TOULOUSE  
FRANCE  
Tel:+33-05-61-11-90-94  
Fax:+33-50-61-48-21-67


### Stable

Designed to match precisely  
your ultrasound equipment's settings.



### Adaptable

Easy to use, and allows  
for a faster biopsy workflow.



### Disposable

Using Leapmed biopsy kits  
reduces greatly the risk of infection.

